

Technical Agreement

www.newarebattery.com



MHW-25-S-16CH Constant Temperature Chamber

(With Integrated Battery Testing System)

Technical Agreement

Neware Technology Limited

Contact: jade@newarebattery.com



Technical Agreement

www.newarebattery.com

Product Name:

Desktop Constant Temperature Chamber

(With Integrated Battery Testing System)



P.S. Image is for reference purpose only.

Model Code	MHW-25-S-16CH		
Application	Constant temperature test of coin cells.		
Prohibitions	Testing or storage of :		
	- flammable, explosive and volatile material samples;		
	- corrosive substances;		
	- strong electromagnetic emission source;		
	- radioactive material samples;		
	- highly toxic substances;		
	- samples that may produce the above substances or objects during		
	testing or storage.		

Dimensions

Nominal Volume	25L
Inner Dimension	W360 mm×D300 mm×H235 mm
Outer Dimension	W440 mm×D580 mm×H410 mm
Net Weight	Around 45kg

Performance

Testing Environment	Room Temperature above 25 $^{\circ}\mathrm{C}$		
	Relative humidity≤85%		
Temperate Range	15∼60℃		
Fluctuation	≤1°C (No load, or during stable temperature)		
Deviation	±2.0°C (No load, or during stable temperature)		
Heating Time	25°C →60°C ≤40 min (No load, average non-linearity)		
Cooling Time	25°C →15°C ≤60 min (No load, average non-linearity)		

Web: www.newarebattery.com Email: jade@newarebattery.com Sales Rep: Jade Zhao



Technical Agreement

www.newarebattery.com

Structure			
	0.1		
	Outer wall material:		
Insulation Envelope	High quality cold rolled steel plate with surface spray paints.		
	Inner wall material: Stainless Steel SUS304		
A	Insulation material: Polyurethane foam		
Air-conditioning Channel	Axial Fan, Semiconductor cooling/heating module		
Standard	Door: insulated tempered glass + Frames		
Configuration	Cell Connection Plate: 2pcs		
Comiguration	Illumination: LED		
Control Panel	Control buttons		
Air Conditioning Unit	Semiconductor cooling/heating module		
Electrical Control	System		
Controller	LED digital display + button controller		
Setting Method	Button controller		
	Forced circulation ventilation.		
Control Method	The system controls the output of the semiconductor		
Control Method	refrigeration/heating module through the PID results, in order to achieve		
	a dynamic balance.		
Communication	Ethernet		
Power Cable			
Cables	Single phase + protective ground wire		
Conditions of Use			
	Good ventilation;		
In stallation Cita	No strong vibration around the device;		
Installation Site	No strong electromagnetic fields around the device;		
	No flammable/explosive/corrosive substances & dust around the device.		
Environmental	Temperature: 25°C±3°C;		
Conditions	Relative humidity: ≤85%;		
COHUILIONS	Atmospheric pressure: 86kPa~106kPa		
Power Supply	Input: AC(220±10%)V /50Hz	Input: AC(110±10%)V /60Hz	
Condition	Max. current: 1A	Max. Current: 2A	
Precautions	Opening the door while testing will cause temperature fluctuations.		

Web: www.newarebattery.com Email: jade@newarebattery.com Sales Rep: Jade Zhao